SAFE GUARD®

24CM PRESSURE ASSISTED DEVICE

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INSTRUCTIONS FOR USE

Read instructions prior to use.

Product not made with natural rubber latex.

DEVICE DESCRIPTION

The SAFE GUARD 24 cm is a single use disposable device. SAFE GUARD has a clear medical grade polyurethane window and bladder, a clear medical grade PVC flexible fill tube, and a pressure sensitive, self-adhesive pelvic backing. A luer valve on the end of the fill tube enables a syringe to be connected to inflate the central bladder with air to provide pressure to the puncture site. The SAFE GUARD pressure assisted device has a sterile dressing with a clear window that facilitates visibility of the access site without removal or manipulation of the device.

INDICATIONS

The indications for use for the SAFE GUARD 24 cm pressure assisted device are to assist in obtaining and maintaining hemostasis. The device is also indicated in the reduction of active compression time in femoral artery cannulation following diagnostic and interventional procedures.

CONTRAINDICATIONS

The adhesive portion of the SAFE GUARD device should not be used over excoriated skin.

CAUTIONS

R Only. Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

CAUTION: With over-inflation i.e., above 40 mL's of air, the bulb may begin to expand radially and could compromise the adhesive properties of the device.

CAUTION: Under-inflation of device could compromise the ability of the device to assist in obtaining and maintaining hemostasis.

REUSE PRECAUTION STATEMENT

For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

PRECAUTIONS

- Use proper aseptic techniques while handling product.
- Do not use if package is damaged
- Inspect device prior to use to verify that no damage has occurred during shipping.

POTENTIAL ADVERSE EFFECTS

Possible adverse effects that may result from the use of this device:

- Hematoma
- Local bleeding
- Arterio-venous fistula or pseudoaneurysm

PRE-HEMOSTASIS, or MANUAL ASSIST TECHNIQUE (MAT)

PLACEMENT OF SAFE GUARD

1. Before adhering SAFE GUARD to the patient, be sure skin is clean and dry. Determine the appropriate angle for SAFE GUARD placement to provide easy access to luer inflate/deflate port and to allow for easy sheath removal.

Note: Placement may require adjustment based on the patient's anatomy, angle of the puncture site, and the presence or absence of a procedural sheath.

2. Consider the point of maximum pulse, anatomy, angle of puncture and direction of flow to determine the appropriate SAFE GUARD position and verify.
3. Pull the procedural sheath back approximately 1" (2.5 cm) so that when SAFE GUARD is adhered to the skin the sheath hub is outside the area of the SAFE GUARD adhesive.

Note: It is recommended that you aspirate the sheath prior to removal to prevent distal embolization from residual clot in sheath.

4. Remove the adhesive backing and place the bulb where you would position your fingers to hold manual compression (for example, in femoral artery procedures, typically the point of maximum femoral pulse). Make sure SAFE GUARD is completely adhered to the skin.
5. Attach and completely engage a standard luer lock syringe to inflate the desired volume a maximum volume of 40 mL's of air into the bulb to apply pressure on the arteriotomy site. Syringe must be completely engaged in the luer to inflate/deflate the bulb. Remove syringe.

Note: Maintain pressure on the plunger while detaching syringe from the SAFE GUARD valve. Observe that the desired pressure is achieved and maintained.






6. Remove sheath, then immediately apply manual compression directly over inflated bulb.
7. Hold manual compression until hemostasis has been achieved.*
- Slowly release manual compression.
- Check distal/proximal pulses to assure flow is maintained.
- Confirm hemostasis by viewing the site through the inflated bulb window.

* Recommendations (MAT only):

Diagnostic patients - minimum 5 minutes

Interventional patients - minimum 10 minutes

8. Per hospital protocol, periodically check the site through the bulb window to confirm hemostasis and to manage the bulb volume and resultant pressure as needed. Continue to check distal/proximal blood flow to assure patency.
9. Deflate bulb every two hours to allow for capillary refill and to assess the site. Re-inflate the bulb if necessary.


- Deflate the bulb by attaching an appropriately sized luer lock syringe to the valve, engage the valve and slowly depress the bulb allowing the syringe to fill with air. Alternatively, remove plunger from the syringe, attach syringe and allow air to slowly release while gently depressing the bulb.

Note: Do not draw negative pressure in the syringe, as this will create a vacuum on the site.

- Prior to discharge of the patient, remove SAFEGUARD and apply sterile dressing per hospital protocol.

POST-HEMOSTASIS TECHNIQUE


- When hemostasis at the access site has been achieved, apply the SAFEGUARD device with the access site visible under the bulb window of the SAFEGUARD device. Consider the point of maximum pulse, anatomy, angle of puncture and direction of flow to determine the appropriate SAFEGUARD position and verify.


Note: Before adhering SAFEGUARD to the patient, be sure skin is clean and dry. Determine the appropriate angle for SAFEGUARD placement to provide easy access to luer inflate/deflate port.

- Attach an appropriately sized standard luer lock syringe to the valve of the SAFEGUARD device.

Note: Syringe must be completely engaged in the luer to inflate/deflate the bulb.

- Inflate the bulb of the SAFEGUARD device with air to the desired volume of air (24cm maximum of 40 mLs) to apply pressure on the arteriotomy site and remove the syringe. Check distal/proximal pulses to assure flow is maintained.

Note: Maintain pressure on the plunger while detaching syringe from the valve of the SAFEGUARD device. Observe that the desired volume is achieved and maintained.

- Per hospital protocol, periodically check the site through the bulb window to assure hemostasis is maintained and the bulb maintains pressure.
- Deflate bulb every two hours and assess the site. Re-inflate the bulb if necessary.
- Deflate the bulb by attaching an appropriately sized luer lock syringe to the valve, engage the valve and slowly depress the bulb allowing the syringe to fill with air. Alternatively, remove plunger from the syringe, attach syringe and allow air to slowly release while gently depressing the bulb.

Note: Do not draw negative pressure in the syringe, as this will create a vacuum on the site.

- Prior to discharge of the patient, remove the SAFEGUARD device and apply sterile dressing per hospital protocol.

Maintain sterile field during application.

Keep away from sunlight

Do not Resterilize

Do not Reuse

Do not use if package is damaged

Catalog Number

Lot Number

Use By

STERILE | R Sterilized Using Gamma

Caution: Consult accompanying document

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Manufacturer: Merit Medical Systems, Inc.
1600 West Merit Parkway, South Jordan, Utah 84095 U.S.A. 1-801-253-1600
U.S.A. Customer Service 1-800-356-3748
Authorized Representative: Merit Medical Ireland Ltd,