Prelude®

Short Sheath Introducer


INSTRUCTIONS FOR USE

PRODUCT DESCRIPTION

The Merit Prelude® Short Sheath Introducer consists of the following components. These components may be packaged in a single pouch or may be packaged separately.

- One (1) Sheath introducer with marker tip and detachable stopcock
- One (1) Vessel dilator
- One (1) Mini guide wire

Some kits could also include an access needle.

- A. Dilator
- B. Detachable stopcock
- C. Side port extension
- D. Removable arterial/venous (red/blue) indicators
- E. Snap clamp
- F. Sheath introducer hub with hemostasis valve
- G. Rotating suture ring
- H. Radiopaque marker tip

INTENDED USE

The Merit Prelude Short Sheath Introducer is intended to provide access and facilitate the percutaneous introduction of various devices into veins and/or arteries while maintaining hemostasis for a variety of diagnostic and therapeutic procedures. This device can also provide access to a native or synthetic graft used for hemodialysis. The side port of the sheath allows adequate flow to perform temporary hemodialysis. This device is not indicated for long term vascular or hemodialysis access.

WARNINGS

- Do not advance the introducer and/or guide wire if resistance is met.
- Do not use a power injector to inject through the side port of the sheath introducer.
- Remove any potential air and/or debris from the device utilizing standard flush protocol prior to attempting infusion through side port of sheath introducer.
- Inspect device prior to use.
- Utilize caution when attempting a second needle/sheath access after the first introducer sheath has been placed to prevent accidental perforation of the first introducer sheath.
- Always flush device appropriately according to hospital protocol to prevent potential clot formation and/or debris accumulation in the fluid path.
- Do not leave the introducer in place for extended periods of time without a catheter or an obturator to support the cannula wall.

CAUTIONS

- Read instructions prior to use
- Store in a cool dry place
- This device is intended for single use only. Do not reuse or resterilize.

POTENTIAL COMPLICATIONS

Potential complications include, but are not limited to:

- Air embolism
- Hematoma
- Perforation or laceration of the vessel wall
- Guide wire embolization
- Risks normally associated with percutaneous diagnostic and/or interventional procedures.
- Inflection
- Bleeding
- Thrombus formation
- Pseudo aneurysm formation
- Vessel spasm

INSTRUCTIONS FOR USE

- Identify the insertion site and prepare the site using proper aseptic technique and local anesthesia as required.
- Remove the Prelude Short Sheath Introducer components from package using proper aseptic technique.
- Flush all components with saline or suitable isotonic solution. The side port of the device has a snap clamp; ensure that clamp is open during flushing.
- After flushing side port, turn stopcock to off position to maintain flush in side port and prevent bleed back upon insertion into the vessel. Close snap clamp on side port extension.
- The side port of the device includes two removable indicators – red indicates arterial and blue indicates venous use. Select the appropriate color and leave in place on the side port tubing; remove and discard the other. If indicators are not necessary, remove and discard both.
- Insert dilator into Prelude Short Sheath Introducer through hemostasis valve and snap into place. Dilator must be securely snapped into place to avoid damage to the vessel.
- Insert appropriate access needle as needed. While holding the access needle, place the flexible end of J end of the guide wire through access needle into the vessel. Note – Refer to product labeling for the guide wire size that is compatible with the system components.

WARNING: Do not advance the guide wire if resistance is met. Determine the cause of resistance before proceeding.

NOTE: If a crossed catheter technique is utilized, the punctures should be greater than 11 cm apart. This distance ensures that during hemodialysis the Prelude Short Sheaths are sufficiently separated for use as arterial and venous conduits.

- Hold guide wire in place and remove access needle. Hold pressure at the site until the introducer/dilator assembly is placed.
- WARNING:** If a needle with a metal cannula is used, do not withdraw the guide wire after it has been inserted because it may damage the guide wire.
- Insert the introducer/dilator assembly over the guide wire into the vessel. Using a rotating motion, advance the introducer/dilator assembly through the tissue into the vessel. Grasp the assembly close to the skin as it is being placed into the vessel to avoid buckling.

- After introducer/dilator assembly has been placed into vessel, detach the dilator from the introducer by bending the dilator hub down slightly (This will un-snap the dilator hub from the introducer cap). While holding the sheath, carefully remove the dilator and guide wire together, leaving the sheath introducer in the vessel.
- Aspirate from the side port extension to remove any potential air. After aspiration, flush the side port with a suitable solution. Stopcock should be turned off to maintain flush in side port. In addition, ensure that clamp is closed to also maintain flush in side port.
- Insert selected device(s) (wires, catheters, etc.) into Prelude Short Sheath Introducer.

NOTE: Hold the sheath in place when inserting, positioning, or removing the devices. Always exchange or remove devices slowly through the sheath.

- Temporarily suture the sheath in place, use the rotating suture ring.
- REMOVAL
 - The sheath may be removed when clinically indicated. Compression on the vessel, above the puncture site, should be started as the sheath is slowly removed. Discard the sheath appropriately.
- NOTE:** Collected fibrin at the tip of the sheath may be aspirated via the side port tubing prior to removal of the sheath.
- Follow standard protocol for achieving hemostasis at the insertion site.
- HEMODIALYSIS
 - If the sheath is to be used for temporary access for hemodialysis, aspirate from the side port extension to remove any potential air. After aspiration, flush the side port with a suitable solution. Stopcock should be turned off to maintain flush in side port. In addition, ensure that clamp is closed to also maintain flush in side port.
 - Patient should be accompanied to dialysis center by appropriate personnel.
 - Following standard protocol, side port of sheath introducer can be connected to hemodialysis circuit.
 - At conclusion of dialysis, refer to step 14 for sheath removal.

Manufacturer:

Merit Medical Systems, Inc. South Jordan, Utah 84095 www.merit.com
U.S.A. 1-801-253-1600 U.S.A. Customer Service 1-800-356-3748

Authorized Representative:

Merit Medical Ireland Ltd, Parkmore Business Park West, Galway, Ireland
European Customer Service by Country: Belgium 0800 72906;
France 0800 916030; Germany 0800 1820871; Ireland 091 7037700;
Neth. 0800 0228184; U.K. 0800 973115

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