


CentrosFLO®

Long-Term Hemodialysis Catheter

- Preshaped curve tips stabilize and center the catheter
- Provides up to **27%** better patency at 3 months than conventional split tip catheter³
- Helps prevent thrombosis, fibrin sheathing & vessel wall occlusion



CentrosFLO®

Long-Term Hemodialysis Catheter

centered on performance

The high performance CentrosFLO® long-term hemodialysis catheter from Merit Medical was designed with performance, safety, and ease of use in mind. The **SELF-CENTERING**, curved-tip catheter has a greater separation between the arterial and venous tips than symmetrical tip catheters. Its design is intended to optimize **CATHETER PATENCY**, and reduce fibrin sheath formation, thrombosis and vessel wall occlusions by keeping the tips of the catheter centered in the vessel and away from the vessel wall. Fibrin sheaths are a frequent cause of catheter malfunction and in a study of removed or exchanged hemodialysis catheters, 76% had fibrin sheaths.^{1, 2, 3} The CentrosFLO was designed for **PERFORMANCE, SAFETY AND EASE OF USE**.

performance

Preshaped Curved Tips

stabilize and center the catheter in the junction of the Superior Vena Cava (SVC) and Right Atrium (RA).


Distal Venous and Arterial Side Holes

reduce intraluminal pressure and the likelihood of the catheter tips sucking up against the vessel wall.

Large Internal Lumen

allows for a flow rate of up to 450 mL/min.

Tunneler with Sleeve



Peelaway Sheath and Dilator

safety

Atraumatic Catheter Tip Design

minimizes puncture or irritation of vessel walls.

Non-Traumatic Tunneler

for more rapid loading of the catheter without risk to catheter tips.

Multiple Catheter Length Configurations

accommodate variations in patient anatomy.

ease of use

Peelaway Sheath and Dilator

insertion for initial placement.

Stiffening Stylet

for over-the-wire placement and exchange.

Polyurethane Material

provides ease of insertion and maintenance.

Stiffening Stylet

CentrosFLO®
Long-Term
Hemodialysis Catheter

study results

In recent studies, the CentrosFLO was compared to two different tunneled dialysis catheters and in both cases the CentrosFLO had statistically better patency. In the first study, the patency of the CentrosFLO was compared to the patency of a symmetric-tip catheter. In the second study, the CentrosFLO's patency was compared to a conventional split-tip catheter when either catheter was exchanged over a wire using the same access site as a replacement for an existing dysfunctional (blood flow rate less than 300mL/min) catheter.

Study A ⁴	Self-Centering, Split-Tip Design (n=38)	Symmetric-Tip Design (n=20)	Statistical Probability
Follow-up Period	3 Months	3 Months	3 Months
3 Month Patency	94.7%	75.0%	P=.04
Mean Clearance (Kt/V)	1.48	1.42	(p=.48)*
Mean Blood Flow (mL/min)	394	382	(p=.27)*

Study B ⁵	Self-Centering, Split Tip Design (n=46)		Conventional Split-Tip Design (n=27)		Statistical Probability	
Follow-up Period	3 Months	6 Months	3 Months	6 Months	3 Months	6 Months
Patency	67.4%	23.9%	40.7%	14.8%	p=.01	p=.01
Mean Blood Flow (mL/Min)	391	388	388	352	(p=.6)*	p=.01


*NDS = Not Statistically Significant, p>0.05 (t-tests)
Note: Mean clearance and mean blood flow did not include data from the patients with no catheter patency and represent values from the follow-up period.

ordering information



complete kit

- 1 15F CentrosFLO Long-Term Hemodialysis Catheter
- 1 Stiffening Stylet
- 1 16F Peelaway Sheath And Dilator
- 1 18G X 2.75" (7 cm) Introducer Needle
- 1 #11 Safety Scalpel
- 1 0.038" X 80 cm J-Tipped Guide Wire
- 2 Adhesive Dressing
- 2 Injection Caps
- 1 12F Vessel Dilator
- 1 14F Vessel Dilator
- 1 Tunneler With Sleeve



catheter only

- 1 15F CentrosFLO Long-Term Hemodialysis Catheter
- 1 Stiffening Stylet
- 2 Injection Caps

All Products Packaged 5 Per Box.

Complete Kit		Quantity Requested		
product code	tip-to-cuff	tip-to-hub	stylet	
CENFP15K	15 cm	20 cm	yes	<input type="text"/>
CENFP17K	17 cm	22 cm	yes	<input type="text"/>
CENFP19K	19 cm	24 cm	yes	<input type="text"/>
CENFP23K	23 cm	28 cm	yes	<input type="text"/>
CENFP27K	27 cm	32 cm	yes	<input type="text"/>
CENFP31K	31 cm	36 cm	yes	<input type="text"/>

Catheter Only		Quantity Requested		
product code	tip-to-cuff	tip-to-hub	stylet	
CENFP15C	15 cm	20 cm	yes	<input type="text"/>
CENFP17C	17 cm	22 cm	yes	<input type="text"/>
CENFP19C	19 cm	24 cm	yes	<input type="text"/>
CENFP23C	23 cm	28 cm	yes	<input type="text"/>
CENFP27C	27 cm	32 cm	yes	<input type="text"/>
CENFP31C	31 cm	36 cm	yes	<input type="text"/>

Ordering Department _____ Customer Signature _____ Date _____

To find out more about the CentrosFLO, call Customer Service at 1-800-35-MERIT.

¹ Crain, MR, Horton, MG, Mewissen, MV. Fibrin sheaths complicating central venous catheters. *AJR*. 1998 Aug;171:341-346.
² Forauer AR, Theoharis GA, Dasika NL. Jugular vein catheter placement: histologic features and development of catheter-related (fibrin) sheaths in a swine model. *Radiol*. 2006 Aug;240(2):427-434.
³ Alomari, Al, Falk, A. The natural history of tunneled hemodialysis catheters removed or exchanged: a single-institution experience. *JVIR*. 2007;18:227-235.
⁴ Balamuthusamy S. Self-centering, split-tip catheter has better patency than symmetric-tip tunneled hemodialysis catheter: single-center retrospective analysis. *Semin Dial*. 2014 Jan 20. Epub 2014 Jan 20.
⁵ Data on file. Exclusions from Study B were patients with catheter-related infections, femoral vein catheters, or exchanges requiring the creation of a new venotomy or exist site. A repeated measures ANOVA was used to assess differences between groups and among times.



Understand. Innovate. Deliver.™

Merit Medical Systems, Inc.
 1600 West Merit Parkway
 South Jordan, Utah 84095
 1.801.253.1600
 1.800.35.MERIT

Merit Medical Europe, Middle East, & Africa (EMEA)
 Amerikalaan 42, 6199 AE Maastricht-Airport
 The Netherlands
 +31 43 358 82 22
 Merit Medical Ireland Ltd.
 Parkmore Business Park West
 Galway, Ireland
 +353 (0) 91 703 733

Merit.com

Austria
 0800 295 374

Belgium
 0800 72 906 (Dutch)
 0800 73 172 (Français)

Denmark
 80 88 00 24

Finland
 0800 770 586

France
 0800 91 60 30

Germany
 0800 182 0871

Ireland (Republic)
 1800 553 163

Italy
 800 897 005

Luxembourg
 8002 25 22

Netherlands
 0800 022 81 84

Norway
 800 11629

Portugal
 308 801 034

Russia
 +7 495 221 89 02

Sweden
 020 792 445

Switzerland
 (Deutsch)
 +41 225180252
 (Français)
 +41 225948000
 (Italiano)
 +41 225180035
 UK
 0800 973 115

401060001/F 2015-12-07