Careflow™ Central Venous Catheter

Careflow™ Arterial Catheter

EN	INSTRUCTIONS FOR USE
FR	MODE D'EMPLOI
IT	ISTRUZIONI PER L'USO
DE	GEBRAUCHSANWEISUNG
ES	INSTRUCCIONES DE USO
PT	INSTRUÇÕES DE UTILIZAÇÃO
NL	GEbruIKSAANWIJZING
SV	BRUKSANVISNING
NO	BRUKSANVISNING
DA	BRUGSANVISNING
EL	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
TR	KULLANMA TALİMATLARI
PL	INSTRUKCJA UŻYTKOWANIA
CS	NÁVOD K POUŽITÍ
BG	ИНСТРУКЦИИ ЗА УПОТРЕБА
HU	HASZNÁLATI ÚTMUTATÓ
RU	ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ
LT	NAUDOJIMO INSTRUKCIJOS
RO	INSTRUCȚIUNI DE UTILIZARE
SK	NÁVOD NA POUŽITIE
HR	UPUTE ZA UPORABU
ET	KASUTUSJUHEND
FI	KÄYTTÖOHJEET
ZH	使用说明

Careflow™ Catheter

INSTRUCTIONS FOR USE

The Merit Careflow™ Catheter is a radio-opaque polyurethane catheter. (Central Venous, Peripheral Venous Catheters and Radial, Femoral Arterial Catheter).

INTENDED USE

The Merit Careflow Central Venous Catheter is an intravenous catheter that is intended for short term use (no longer than 30 days) to access the human circulatory system via insertion through either the Internal Jugular or Subclavian Vein using Seldinger technique, where the catheter tip resides in the superior vena cava. The Central Venous Catheter is intended for infusion of drugs, Total Parenteral Nutrition (TPN) fluids, large volume infusions, repeated blood sampling and to monitor the central venous pressure.

The Merit Careflow Arterial Catheter is intended for short term use (no longer than 30 days) to assess the radial and femoral arteries for repeated blood sampling and to allow invasive blood pressure measurement.

Note: Use only as directed by a physician.

A full list of components is shown on the package lid.

GENERAL WARNINGS

- Physicians must be familiar with the complications associated with central venous catheterisation, i.e. vessel perforation, air embolism, catheter embolism, pleural and mediastinal damage, septicemia, thrombosis and cardiac tamponade secondary to vessel wall or atrial perforation.
- Complications are associated with right atrial and inadvertent right ventricular catheterisation. Physicians must be aware of these complications before advancing the catheter beyond the depth required for normal vena cava placement. Do not advance the catheter past this depth unless procedure requires right atrial placement. If catheter is advanced beyond normal vena cava placement depth, monitor electrocardiogram during insertion and confirm final position by chest X-Ray.
- It is recommended that patients be placed in a slight Trendelenburg position during insertion procedure to reduce the risk of air embolism.
- We recommend that Careflow™ Catheter lumens (Central Venous, Peripheral Venous Catheters and Radial, Femoral Artery Catheter) are flushed with sterile saline solution prior to catheter insertion.
- All catheter placements must be inspected for flow rates, security of dressing and security of luer connections.
- To reduce or eliminate the potential for catheter migrations, we advise that every catheter placement be secured by suturing at the eyelets of the junction boot housing/hub and that, where the use of the secondary fixation device is necessary, it should be used as additional support and not the only means of fixation. Additionally, the security of catheter fixation and position of the catheter tip should be checked throughout use.
- When removing dressings at or close to catheter sites, care must be taken to avoid severing the catheter.
- Acetone must not come into contact with the catheter as the material may weaken and this may result in leakage or aspiration.
- Exposure of product componentry to topicals containing alcohol is not recommended.
- Do not attempt to re-insert a partially or completely withdrawn introducer cannula.
- Percutaneous puncture of a central vein may be contraindicated in patients with pulmonary hypertension.
- Use of a syringe smaller than 5 mL to irrigate or de-clot an occluded catheter may cause intraluminal or catheter rupture.
- Syringes are supplied for blood aspiration only.
- Luer connections: as standard practice the security of luer connections must be checked routinely.
- Physicians should be aware that central venous and arterial catheters are intended for use no longer than 30 days.

WARNINGS - SELDINGER TECHNIQUE

- Do not withdraw the guidewire against needle bevel, as this increases the risk of severing the guidewire.
- During insertion do not reinsert a partially or completely withdrawn needle into the cannula.
- Ensure the flexible end of the guidewire is advanced into the vein.
- Ensure the guidewire moves freely in the needle introducer.
- A firm grip must be maintained on the guidewire at all times.
- When using the 'J' wire straightener maintain a firm grip on the plastic sleeve.
- Ensure the dilator is removed prior to catheter advancement.
- The moveable suture devices are designed as additional support and must not be used as the only means of fixation.
- Potential for guidewire breakage. Although the incidence of guidewire breakage is extremely uncommon, physicians must be aware of the potential of guidewire breakage if undue force is applied to the wire. If resistance is met when attempting to remove the guidewire after central venous placement, the wire may be kinked within the area of the catheter tip and the vessel. Undue force may cause the wire to break. If resistance is encountered, withdraw the catheter relative to the guidewire (2-3 cm) and attempt to remove the wire. If resistance is still apparent remove the wire and the catheter simultaneously.
- Physicians should be aware that the guidewire can pick up material from the vein. This may prevent the guidewire from being withdrawn through the catheter.
- Do not force the guidewire. If resistance is met, carefully withdraw the guidewire and re-attempt insertion.
- If femoral approach is used, place patient in supine position for insertion procedure.

SUGGESTED PROCEDURE SELDINGER TECHNIQUE (VENOUS ACCESS)


- Prepare the insertion site using the full aseptic technique required for venous access.
- The guidewire may be introduced either by using a thin walled steel needle or introducer cannula.
- Feed the desired flexible end of the guidewire through the introducer into the vein. If a 'J' wire is to be used, the 'J' can be straightened prior to insertion by using the plastic insertion sleeve.
- Remove the introducer.
- A vessel dilator may be used to enlarge the cutaneous puncture site. If further enlargement of puncture site is necessary, use a scalpel.
- When using multi lumen catheters, lumens other than distal should be flushed and attached to desired fluid administration set or, alternatively, heparin locked as standard hospital practice.
- Pass the distal tip of the catheter over the guidewire. (Sufficient guidewire length must remain exposed at the hub end of the catheter to maintain a firm grip on the guidewire.) Grasping the

catheter near the skin, push the catheter into the vein with a slight twisting motion and advance into the final position.

- Hold the catheter in position, withdraw the guidewire and aspirate with a syringe to ensure correct placement.
- The distal hub should be connected to the appropriate fluid administration set. If a hub with Floswitch™ is being used this can now be switched off.
- The catheter may now be secured by suturing eyelets of the junction boot housing/hub to the skin.
- Apply sterile dressing as appropriate.

SECONDARY FIXATION DEVICE (5Fr-8.5Fr CATHETERS)


Note: This device is found in certain models only.

- Position the secondary fixation device on catheter. To close the device, press down on each wing. A 'click' sound confirms that the device is secured and cannot be moved easily.
- To ensure fixation, suture at each eyelet on the wings of the Secondary Fixation Device.
- To open the device, press down on the 2nd hinge, pull up on each wing.

SECONDARY FIXATION DEVICE (7Fr AND 9.5Fr CATHETERS)


Note: This device is found in certain models only.

- Ensure the Secondary Fixation Device is in place at the point of suture.
- To close the device, snap the top plate in place.
- To ensure fixation, suture at each eyelet on the wings of the Secondary Fixation Device.
- To open the device, hold the lock device down and pull the plate up, using the lever.

VENAGUIDE™

Note: This device is found in certain models only.

- Release guidewire by removing the guidewire cap.
- Straighten guidewire 'J' by retracting into introducer system with thumb.
- Insert into hub of introducer needle and advance guidewire into vein. Advance to required depth. **WARNING:** The included guidewire does not have depth markings. Electrocardiogram (ECG) monitoring or Ultrasound guidance, or a combination of both is recommended to prevent under insertion/ over insertion of the guidewire.

GUIDING SYRINGE

Note: This accessory is found in certain models only.

- This device is used for introducing guidewires for use with Merit central or peripheral venous catheters.
- Insert introducer needle attached to Guiding Syringe into vessel and aspirate.
 - Feed the desired flexible end of the guidewire through the hole in rear of Guiding Syringe plunger into the vein. If a 'J' wire is to be used, straighten guidewire 'J' by retracting into introducer system with thumb.
 - Hold the guidewire in place and remove introducer needle and Guiding Syringe. **WARNING:** Do not aspirate with guidewire in place or air may enter the syringe. **Caution:** Do not reinfuse blood to minimize the risk of blood leakage from rear of syringe.

FLOSWITCH™ INTRODUCER

This device is used for introducing guidewires for use with Merit central or peripheral venous catheters.


Using a Floswitch™ introducer cannula, puncture the vein, aspirate and advance the cannula into the vein and remove the needle. At this stage the Floswitch™ may be switched off. The switch must be opened again to facilitate guidewire or catheter passage.

FLOSWITCH™ LUER LOCK ATTACHMENT

This device may be attached to a luer hub of a central or peripheral venous catheter.

- As standard practice the catheter should be secured to the skin.
- To attach the Floswitch™, insert the Floswitch™ luer into the catheter hub. Tighten the rotating collar ensuring that 'hand tight' connection has been made and that the catheter hub is fully advanced into the Floswitch™ collar (see figure 1.1).
- The Floswitch™ should be fixed in position by either suturing or taping the wings. The recommended procedure for taping is shown in figure 2.
- The black markings indicate flow status. When visible the catheter is open and there is free flow. When covered, the catheter is closed.


Figure 1. Floswitch™ Inspection Procedure

1.1 Catheter hub fully advanced into Floswitch collar

1.2 Inspection of connection necessary

Figure 2. Floswitch Taping Procedure


WARNINGS - FLOSWITCH™ PRODUCTS

1. The Floswitch™ should not be switched off (black marks covered) before the needle has been fully withdrawn. This applies to Floswitch™ introducers only.
2. The Floswitch™ must not be switched off (black marks covered) until the guidewire has been fully withdrawn.
3. When the catheter is not in use for infusion or aspiration, the Floswitch™ must be switched off and a suitable luer cap locked into the hub.
4. When using an intermittent injection bung attached to a Floswitch™, only use short needles. Ensure Floswitch™ is in 'ON' position prior to injecting or aspirating. Do not insert the needle more than 8 mm into the Floswitch™.
5. As standard practice, the security of the luer connection must be checked routinely. This is essential when lubricious substances such as Intralipids are being used. This applies to Floswitch™ luer lock attachment only.

FUTURA® SAFETY SCALPEL

Note: This device is found in certain models only.

1. Push forward with thumb until blade locks.

2. Retract blade before passing scalpel or disposing scalpel. To retract apply pressure here.

3. Scalpel is now ready for passing or disposal.

WARNING: Used sharps are contaminated. Handle carefully. Dispose scalpel in accordance with all applicable local and national laws and regulations.

STORAGE CONDITIONS

Store in cool dry place away from direct sunlight.

STERILE AND NON-PYROGENIC in unopened, undamaged package. For single use only. Check integrity of the individual package before use. Dispose of product after use. Do not resterilize.

Re-use may lead to infection or other illness/ injury.

For reordering information or assistance please contact local representative.

	Do Not Use If Package is Damaged
Rx ONLY	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
	Do not Resterilize
	Single Use
	Does not Contain DEHP, DIBP, DBP, BBP
	Non-Pyrogenic
	Sterilized using Ethylene Oxide
	Caution: Consult accompanying documents.
	Use by date: YYYY-MM-DD
	Authorized Representative in European Community
	Manufacturer
	Catalog number
	Batch code

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Merit Medical Singapore Pte. Ltd.,
Merit Medical Ireland Ltd.