

# MERIT MARQUIS<sup>®</sup> SERIES STOPCOCKS



SHERLOCK<sup>®</sup>  
CONNECTORS

AIRLESS ROTATORS

CLEAR INNER LUMEN

SHARKS FIN HANDLE





# Merit Medical

## Marquis® Series Stopcocks



S1SNC — Standard pressure, one-way with swivel nut connector

M1SNC — Medium pressure, one-way with swivel nut connector

H1FLC — High pressure, one-way with fixed luer connector

H1RC — High pressure, one-way with rotating connector

S3SNC — Standard pressure, three-way with swivel nut connector

S4SNC — Standard pressure, four-way with swivel nut connector

M3SNC — Medium pressure, three-way with swivel nut connector

H3SNC — High pressure, three-way with swivel nut connector

M3RRC — Medium pressure, three-way with right rotating connector

H3RRC — High pressure, three-way with right rotating connector

Other configurations are available, ask your customer service representative.

S3FLC — Standard pressure, three-way with fixed luer connector

S4FLC — Standard pressure, four-way with fixed luer connector

M3FLC — Medium pressure, three-way with fixed luer connector

H3FLC — High pressure, three-way with fixed luer connector

L3SNC — Large bore, low pressure, three-way with swivel nut connector

L4SNC — Large bore, low pressure, four-way with swivel nut connector

*The Marquis® Series Stopcock is another innovative accessory product from Merit. Our larger, easy-grip handle makes a Marquis Stopcock effortless to grip and turn. Our patented airless rotator prevents air bubbles while the Sherlock® Connector makes disconnection noticeably easier.*

*Using the same state-of-the-art design that made Merit's Manifold the industry standard, the Marquis Series Stopcock offers unsurpassed feel and ease-of-use. Available in pressure ratings of 200, 500 and 1050 PSI with the easiest turning handles on the market. Of course, the best test of any Marquis Stopcock is a direct comparison with any competing product. Once you've tried them all, you'll know why Marquis Series Stopcocks are truly "a cut above the rest."*

### TECHNICAL SPECIFICATIONS

Body Material:	Polycarbonate
Handle Material:	Delrin
Pressure Rating:	
Large bore:	50 PSI (3.45 BAR)
Standard:	200 PSI (13.8 BAR)
Medium:	500 PSI (34.5 BAR)
High:	1050 PSI (72.45 BAR)
Primary Handle I.D. (Large bore)	.079 inches (2.0 mm)
	.120 inches (2.5mm)
Packaging:	50 units per box, 4 boxes per case



Merit Medical Systems, Inc. ■ 1600 West Merit Parkway ■ South Jordan, Utah 84095 ■ 1-801-253-1600 ■ 1-800-35-MERIT  
Merit Medical EUROPE, AFRICA AND THE MIDDLE EAST (EMEA) ■ Amerikalaan 42, 6199 AE Maastricht-Airport ■ The Netherlands ■ Tel: +31 43 358 82 22  
Merit Medical Ireland Ltd. ■ Parkmore Business Park West ■ Galway, Ireland ■ +353 (0) 91 703 733  
Free phone for specific country: Austria 0800 295 374 ■ Belgium 0800 72 906 (Dutch) 0800 73 172 (French) ■ Denmark 80 88 00 24 ■ France 0800 91 60 30  
Finland 0800 770 586 ■ Germany 0800 182 0871 ■ Ireland (republic) 1800 553 163 ■ Italy 800 897 005 ■ Luxembourg 8002 25 22 ■ Netherlands 0800 022 81 84  
Norway 800 11629 ■ Sweden 020 792 445 ■ UK 0800 973 115

www.merit.com

400177018/D ID 102612