

Bearing nsPVA[®]

Embolization Particles

CATHETER COMPATIBILITY CHART



Merit Maestro[®]
Microcatheter

Merit Pursue[™]
Microcatheter

SwiftNINJA[®]
Steerable Microcatheter



PART NUMBER		V100EP	V200EP	V300EP	V400EP	V600EP	V800EP	V1100EP
PARTICLE SIZE		45-150 µm	150-250 µm	250-355 µm	355-500 µm	500-710 µm	710-1000 µm	1000-1180 µm
CATHETER	INNER DIAMETER							
Angiographic Catheter	0.046" (1.17 mm)							
Angiographic Catheter	0.040" (1.02 mm)							
EmboCath [®] Plus Microcatheter	0.028" (0.71 mm)							
Merit Maestro [®] Microcatheter 2.9F	0.027" (0.69 mm)							
Merit Maestro Microcatheter 2.8F	0.024" (0.62 mm)							
SwiftNINJA [®] Steerable Microcatheter 2.4F	0.021" (0.54 mm)							
Merit Maestro Microcatheter 2.4F	0.020" (0.52 mm)							
Merit Pursue [™] Microcatheter 2.0F	0.020" (0.50 mm)							
Merit Maestro Microcatheter 2.1F	0.018" (0.46 mm)							
Microcatheter	0.017" (0.43 mm)							
Merit Pursue Microcatheter 1.7F	0.016" (0.40 mm)							

Above catheter compatibility information based on internal testing.

Before using refer to Instructions for Use for indications, contraindications, warnings, precautions, and directions for use.



Understand. Innovate. Deliver.™

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1.801.253.1600
1.800.35.MERIT

Merit Medical Europe, Middle East, & Africa (EMEA)
Amerikalaan 42, 6199 AE Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733

merit.com

