



*Spend a little...  
Save a lot.*

*CT transfer sets are used to transfer contrast  
or saline into a sterile injector syringe.*

## **CT TRANSFER SETS**

## CT TRANSFER SETS

*CT transfer sets provide a safe and effective way to decrease turn around time in busy procedure rooms. By utilizing larger contrast bottles, you not only save money on contrast, but also precious time during set-up and post procedure cleanup.*

### FEATURES & BENEFITS



#### CONTRAST TRANSFER SET

- Large lumen spike with auto vent provides high flow
- Wide base unit on piercing spike adds insertion support
- Low priming volume
- Straight through design provides unimpeded flow
- Interchangeable – can be used with either contrast or saline



#### SALINE TRANSFER SET

- Longer, sharper, spike tip allows easier insertion into I.V. solution bags
- Blue stripe tubing is a solution identification reminder
- Interchangeable – can be used with either contrast or saline




#### LARGE BORE CONTRAST TRANSFER SET

- Large bore piercing spike
- 20" large bore tubing
- Large wing swabbable valve



#### SWABBABLE VALVE

- Valve is automatically in the closed position – opens when a syringe tip is inserted into the valve, and closes when you remove syringe tip
- Large wings provide added grip during attachment to syringe
- A quick swab to the valve is all you need to load the next syringe which equals less time spent getting ready for procedures
- Easy access with either slip luer or luer lock syringes



### ORDERING INFORMATION

Catalog Number	Description
<b>TS120</b>	CT Contrast Transfer Set with filtered Sosa spike
<b>TS220</b>	CT Saline Transfer Set with blue stripe tubing and longer spike
<b>TS520L</b>	CT Large Bore Contrast Transfer Set

Can be used with single or dual chamber injectors.  
Packaged: Sterile Single Pack - 50 units/box



Merit Medical Systems, Inc. • 1600 West Merit Parkway, • South Jordan, Utah 84095 • 1-801-253-1600 • 1-800-35-MERIT

Merit Medical Europe • Horsterweg 24 • 6199 AC MAASTRICHT- AIRPORT • The Netherlands • +31 43 35 88 222

Merit Medical Ireland Ltd. • Parkmore Business Park West • Galway, Ireland • +353 (0) 91 703 733

Free phone from specific country:

Belgium 0800 72906 • France 0800 91 60 30 • Germany 0800 182 0871 • Luxembourg 0800 22522 • Netherlands 0800 022 8184 • UK 0800 973 115

www.merit.com

400873003/A ID 090707