



MERIT - MINI ACCESS KITS
MAK - VASCULAR

MERIT MAK™ - MINI ACCESS KITS

The Merit MAK kit provides access to the vascular system and facilitates the placement of 0.035" (0.89 mm) or 0.038" (0.97 mm) guide wires. The MAK kits allow you to start small, minimising bleeding and trauma to the surrounding tissue. Kits are available in 4F and 5F and include the coaxial introducer/dilator pair, 21 gauge needle and 0.018" (0.46 mm) guide wire. For your convenience, all kit components are available separately. Merit MAK guide wires are available in many configurations to meet specific clinical needs. Guide wires are packaged in a fully enclosed hoop to ensure they are protected and ready for use.

KIT INCLUDES:

- Coaxial Introducer/Dilator Pair
- 21G Needle
- Guide Wire



Catalogue Number	French Size	Length	Guide Wire Diameter/Length	Guide Wire Core/Tip	Needle Gauge/Length	Qty/Box
------------------	-------------	--------	----------------------------	---------------------	---------------------	---------

STANDARD 7 CM NEEDLE KITS

MAK401	4	10 cm	0.018" (0.46 mm) x 40 cm	SS/SS	21G x 7 cm	10
MAK501	5	10 cm	0.018" (0.46 mm) x 40 cm	SS/SS	21G x 7 cm	10
MAK401N	4	10 cm	0.018" (0.46 mm) x 40 cm	NT/PT	21G x 7 cm	10
MAK401NPD	4	10 cm	0.018" (0.46 mm) x 40 cm	NT/PD	21G x 7 cm	10
MAK501N	5	10 cm	0.018" (0.46 mm) x 40 cm	NT/PT	21G x 7 cm	10
MAK501NPD	5	10 cm	0.018" (0.46 mm) x 40 cm	NT/PD	21G x 7 cm	10
MAK501N-15	5	15 cm	0.018" (0.46 mm) x 60 cm	NT/PT	21G x 7 cm	10
MAK405	4	10 cm	0.018" (0.46 mm) x 40 cm	SS/PT	21G x 7 cm	10
MAK405PD	4	10 cm	0.018" (0.46 mm) x 40 cm	SS/PD	21G x 7 cm	10
MAK505	5	10 cm	0.018" (0.46 mm) x 40 cm	SS/PT	21G x 7 cm	10
MAK505PD	5	10 cm	0.018" (0.46 mm) x 40 cm	SS/PD	21G x 7 cm	10

ECHO-ENHANCED 7 CM NEEDLE KITS

MAK401-E	4	10 cm	0.018" (0.46 mm) x 40 cm	SS/SS	21G x 7 cm	10
MAK501-E	5	10 cm	0.018" (0.46 mm) x 40 cm	SS/SS	21G x 7 cm	10
MAK401N-E	4	10 cm	0.018" (0.46 mm) x 40 cm	NT/PT	21G x 7 cm	10
MAK401NPD-E	4	10 cm	0.018" (0.46 mm) x 40 cm	NT/PD	21G x 7 cm	10
MAK501N-E	5	10 cm	0.018" (0.46 mm) x 40 cm	NT/PT	21G x 7 cm	10
MAK501NPD-E	5	10 cm	0.018" (0.46 mm) x 40 cm	NT/PD	21G x 7 cm	10
MAK405-E	4	10 cm	0.018" (0.46 mm) x 40 cm	SS/PT	21G x 7 cm	10
MAK505-E	5	10 cm	0.018" (0.46 mm) x 40 cm	SS/PT	21G x 7 cm	10

ECHO-ENHANCED 4 CM NEEDLE KITS

MAK401-PED	4	10 cm	0.018" (0.46 mm) x 40 cm	SS/SS	21G x 4 cm	10
MAK501-PED	5	10 cm	0.018" (0.46 mm) x 40 cm	SS/SS	21G x 4 cm	10
MAK401N-PED	4	10 cm	0.018" (0.46 mm) x 40 cm	NT/PT	21G x 4 cm	10
MAK405-PED	4	10 cm	0.018" (0.46 mm) x 40 cm	SS/PT	21G x 4 cm	10
MAK501N-PED	5	10 cm	0.018" (0.46 mm) x 40 cm	NT/PT	21G x 4 cm	10
MAK505-PED	5	10 cm	0.018" (0.46 mm) x 40 cm	SS/PT	21G x 4 cm	10

SS = Stainless Steel PT = Platinum Tip NT = Nitinol PD = Palladium Tip



Understand. Innovate. Deliver.™



European Corporate Headquarters, Maastricht, The Netherlands

Merit.com

Distribuido por:



ventas@cardioparts.com • www.cardioparts.com
Tels. (55) 9183-3933 (81) 8357-0779 (81) 8357-1903



Merit Medical Europe, Middle East & Africa

Amerikalaan 42, 6199 AE Maastricht-Airport, Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.

Parkmore Business Park West, Galway, Ireland
+35 39 170 37 33

Merit Medical France

136 Ave. Joseph Kessel, 78960 Voisins Le Bretonneux, France
+33 (0)1 30 57 24 24

Merit Medical Systems, Inc.

1600 West Merit Parkway, South Jordan, UT, 84095 USA
+1 801 253 1600

Asia, Pacific Rim, Latin America, Caribbean & Canada

+1 801 208 4438

Austria

0800 295 374

Belgium

0800 72 906 (Dutch)
0800 73 172 (Français)

Denmark

80 88 00 24

Finland

0800 770 586

France

0800 91 60 30

Germany

0800 182 0871

Ireland (Republic)

1800 553 163

Italy

800 897 005

Luxembourg

8002 25 22

Netherlands

0800 022 81 84

Norway

800 11629

Portugal

308 801 034

Russia

+7 495 221 89 02

Sweden

020 792 445

Switzerland

+41 225180252
(Deutsch)

+41 225948000
(Français)

+41 225180035
(Italiano)

UK

0800 973 115