


MERIT MANIFOLDS[®]

When You Need Confidence and Control

246 STYLES AND PRESSURE RATINGS



MERIT[®] MANIFOLD

...with numerous options at your fingertips.
The Merit Manifold gives you options, confidence, and control.



DEVOS MANIFOLD™
Integral Check Relief Valve Feature - Now... self filling syringe and directed flow control with our new manifold with integral check relief valve.

CONTROL DEBUBBLING
with clear, polycarbonate bodies.




CONTROL VOLUME
with large bore (.093"/.25cm) inner lumen throughout, including rotator.

CONTROL ORDERING OPTIONS
with 246 styles and pressure ratings of both 200 and 500 psi.

CONTROL OPERATION
with 6" (15.6cm) extended rotators.



CONTROL ENVIRONMENT
with leak free, smooth turning, low torque handles.



The first 6 alphanumeric characters are always present.
Manifolds with integral check relief valve are only available in 3-port, half body, off handle, right orientation.

*Five port manifolds are only available in BLOCK BODY – ALL FEMALE PORTS.
**When referring to port orientation, remember that LEFT and RIGHT indicates that there is a rotator present. ALL FEMALE implies that there is no rotator.

