

*It's all in
the TOUCH*



basix TOUCH™


The latest innovation from
the world leader in inflation
device technology

basixTOUCH™

Since 1987, Merit Medical has been designing and improving inflation devices to meet the need for faster inflation, greater capacity, versatility, and reliability.

FAST

The average procedure may require multiple inflations. With deflation speed critical, no one wants to spend more time than necessary turning handles or prepping syringes. The basixTOUCH delivers the speed you need to expand the capacity of your lab.



Quick-release handle for easy release even from high-pressure inflations, accelerating deflation time

Ergonomic handle and new threads ensure that high pressures can be reached rapidly – up to 35 atm with less than 3 rotations of the handle

One-handed preparation





One-handed preparation and priming can reduce set up time

VERSATILE

Today's clinicians have to be prepared to deliver a wide variety of stents and balloons. The versatility of the basixTOUCH can provide you the extra performance you need.

30 mL syringe barrel accommodates additional fluid for large balloons.

35 atm pressure capacity allows inflation of high pressure interventional balloons.



ROBUST AND RELIABLE

Busy laboratories cannot afford to have an inflation device fail. With stronger threads, nut, and handle, the basixTOUCH provides expanded capability with the security and confidence in performance that comes from using a Merit product.

The basixTOUCH – Expanding your capacity with the innovation, quality, and reliability that Merit customers have enjoyed for over 25 years.

basixTOUCH™ *It's all in the TOUCH*

ORDERING INFORMATION

Catalog Number

IN8130 basixTOUCH inflation device

Packaged 5 per box.

Order the basixTOUCH with your choice of Merit Medical's market-leading line of hemostasis valves and kit configurations:

IN8112 with Access-9™ hemostasis valve, insertion tool and torque device

IN8152 with AccessPLUS™ hemostasis valve, insertion tool and torque device

IN8302 with Honor® hemostasis valve, insertion tool and torque device

IN8403 with MBA™ hemostasis valve, insertion tool and torque device



Merit Medical Systems, Inc. • 1600 West Merit Parkway • South Jordan, Utah 84095 • 1-801-253-1600 • 1-800-35-MERIT www.merit.com
Merit Medical EUROPE, AFRICA AND THE MIDDLE EAST (EMEA) Amerikalaan 42, 6199 AE Maastricht-Airport, The Netherlands Tel: +31 43 358 82 22
Merit Medical Ireland Ltd. • Parkmore Business Park West, Galway, Ireland • +353 (0) 91 703 700
Free phone for specific country: Austria 0800 295 374 • Belgium 0800 72 906 (Dutch) 0800 73 172 (French) • Denmark 80 88 00 24
France 0800 91 60 30 • Germany 0800 182 0871 • Ireland (republic) 1800 553 163 • Luxembourg 8002 25 22
Netherlands 0800 022 81 84 • Sweden 020 792 445 • UK 0800 973 115

U.S. and other foreign patents pending.

402843001/A ID 011013